

Denudation and Manipulation SYNGA pipettes for IVF

CE certified, sterile, MEA tested, non toxic, safe for living cells.

CATALOGUE

I. Classic System

SG pipettes without bulbs (bulbs sold separately):

SG Pipette, sterile, pack: 5x10 pipette - 50 pipettes in a pack

Denudation

Short – functional part 14 mm

SG-130S-50	inner Ø 130µm, red
SG-140S-50	inner Ø 140µm, green
SG-150S-50	inner Ø 150µm, white
SG-160S-50	inner Ø 160µm, white

Denudation

Long – functional part 29 mm

SG-130L-50	inner Ø 130µm, red
SG-140L-50	inner Ø 140µm, green
SG-150L-50	inner Ø 150µm, white
SG-160L-50	inner Ø 160µm, white



Manipulation and Blastocyst transfer

SG-180-50	inner Ø 180µm, yellow
SG-250-50	inner Ø 250µm, blue

Classic Bulbs

Bulbs - sterile, 4x5 pcs, sealed in a blister pack by 5 pieces, color may be used to indicates the diameter

A. XtraSoft

SG-012-20XG	green
SG-012-20XY	yellow
SG-012-20XB	blue
SG-012-20XR	red
SG-012-20XW	white

B. Soft

SG-012-20SG	green
SG-012-20SY	yellow
SG-012-20SB	blue
SG-012-20SR	red
SG-012-20SW	white

VI Bulbs

Bulbs with reduced inner volume - sterile, 4x5 pcs, sealed in a blister pack by 5 pieces

SG-VI-20XG	green
SG-VI-20XY	yellow
SG-VI-20XB	blue
SG-VI-20XR	red
SG-VI-20XW	white




SG-021-20SG



SG-VI-20XY

SYNGA s.r.o.

Voděradská 2552, 251 01 Říčany u Prahy, Czech republic
mail@synga.cz, www.synga.cz



SYNGA
TOUCHING LIFE GENTLY

Accessories

SG-007-CUT-R	SG Cutter. Storage of unopened sealed pipettes and the cutter tool to be attached to the right side wall of the IVF cabinet
SG-007-CUT-L	SG Cutter. Storage of unopened sealed pipettes and the cutter tool to be attached to the left side wall of the IVF cabinet
SG-011-1	SG Rack. For storage of up to 6 open pipettes with bulbs in the working area in horizontal or vertical position
SG-021-3	SG Pen Holder VI
SG-021-4	SG Pen Holder C



SG-007-CUT-R



SG-011-1



SG-021-4



SG-021-3

SYNGA s.r.o.

Voděradská 2552, 251 01 Říčany u Prahy, Czech republic
mail@synga.cz, www.synga.cz

II. Premounted system - bulb&pipette

SG Pipette pack: 4x5 pipette with bulb attached, blister of 5 pieces, 4 in a box, sterile individually packed

A. Classic Bulb - short

SG-130S-20X-BL	inner Ø 130µm, red bulb, short pipette
SG-130L-20X-BL	inner Ø 130µm, red bulb, long pipette
SG-140S-20X-BL	inner Ø 140µm, green bulb, short pipette
SG-140L-20X-BL	inner Ø 140µm, green bulb, long pipette
SG-150S-20X-BL	inner Ø 150µm, white bulb, short pipette
SG-150L-20X-BL	inner Ø 150µm, white bulb, long pipette
SG-180-20X-BL	inner Ø 180µm, yellow bulb, short pipette
SG-250-20X-BL	inner Ø 250µm, blue bulb, short pipette

B. VI Bulb

SG-130S-20XVI-BL	inner Ø 130µm, red bulb, short pipette
SG-130L-20XVI-BL	inner Ø 130µm, red bulb, long pipette
SG-140S-20XVI-BL	inner Ø 140µm, green bulb, short pipette
SG-140L-20XVI-BL	inner Ø 140µm, green bulb, long pipette
SG-150S-20XVI-BL	inner Ø 150µm, white bulb, short pipette
SG-150L-20XVI-BL	inner Ø 150µm, white bulb, long pipette
SG-180-20XVI-BL	inner Ø 180µm, yellow bulb, short pipette
SG-250-20XVI-BL	inner Ø 250µm, blue bulb, short pipette



Accessories

SG-31-1	SG Magazine. Storage of premounted pipettes with bulbs, to be attached on the wall
---------	---

SYNGA s.r.o.

Voděradská 2552, 251 01 Říčany u Prahy, Czech republic
mail@synga.cz, www.synga.cz

III. SAS - Synga Adhering System

Classic Bulbs with Adhering ball

Classic Bulbs including metal self Adhering ball – in blister pack by 5 pieces:

A. XtraSoft

SA-012-5XG	green
SA-012-5XY	yellow
SA-012-5XB	blue
SA-012-5XR	red
SA-012-5XW	white

B. Soft

SA-012-5SG	green
SA-012-5SY	yellow
SA-012-5SB	blue
SA-012-5SR	red
SA-012-5SW	white



Accessories


SA-011-1	SA Magnetic Rack , for storage of up to 6 open pipettes with SAS bulbs in the working area in horizontal or vertical position.
SA-021-2	SA Pen Holder C , Conical with magnetic dot - For manual holding with hole for push pull bulb and convenient manipulation, One sided, Dural, used with classical bulbs without metal ball
SA-008-R	SA M dots Wall Rack for Adhesive hanging up to 3 open pipettes with SAS bulbs or with a SAS Pen Holder. - Right side
SA-008-L	SA M dots Wall Rack for Adhesive hanging up to 3 open pipettes with SAS bulbs or with a SAS Pen Holder. - Left side
SA-003	SA Base M dot Adhesive hang Rack for storage of up to 3 open pipettes with SAS bulbs or pen holders



SA-008-L



SA-003



SYNGA
TOUCHING LIFE GENTLY

IV. Synga Flex

Plastic pipettes, sterile, for use with Stripper, 5x10 pipettes, OD = 0,9mm

SF-140-50	inner Ø 140µm
SF-180-50	inner Ø 180µm
SF-250-50	inner Ø 250µm
SF-600-50	inner Ø 600µm



SYNGA s.r.o.

Voděradská 2552, 251 01 Říčany u Prahy, Czech republic
mail@synga.cz, www.synga.cz

V. Samples

SG Pipettes and Bulbs, 10x pipette + 1x Bulb S + 1x Bulb XS + 1x Bulb VI

SG-130S-10XS+VI	inner Ø 130µm, red bulbs, short pipette
SG-140S-10XS+VI	inner Ø 140µm, green bulbs, short pipette
SG-140L-10XS+VI	inner Ø 140µm, green bulbs, long pipette
SG-150S-10XS+VI	inner Ø 150µm, white bulbs, short pipette
SG-150L-10XS+VI	inner Ø 150µm, white bulbs, long pipette
SG-160S-10XS+VI	inner Ø 160µm, white bulbs, short pipette
SG-180-10XS+VI	inner Ø 180µm, yellow bulbs, short pipette
SG-250-10XS+VI	inner Ø 250µm, blue bulbs, short pipette



Premounted – pipettes with bulb attached, blister of 5 pieces

SG-130S-5XVI	inner Ø 130µm, red bulb, short pipette
SG-130L-5XVI	inner Ø 130µm, red bulb, long pipette
SG-140S-5XVI	inner Ø 140µm, green bulb, short pipette
SG-140L-5XVI	inner Ø 140µm, green bulb, long pipette
SG-150S-5XVI	inner Ø 150µm, white bulb, short pipette
SG-150L-5XVI	inner Ø 150µm, white bulb, long pipette
SG-180-5XVI	inner Ø 180µm, yellow bulb, short pipette
SG-250-5XVI	Inner Ø 250µm, blue bulb, short pipette