



TPC Micropipettes

- Micropipettes since 2000
- No compromise on quality
- Designed for optimal outcome and minimal oocyte degeneration



ICSI Micropipettes

| Type | Reference no. | Available angles | Distance tip to bend | Bevel | Inner diameter | Other |
|-----------------------|---------------|------------------|----------------------|-------|----------------|---------------|
| Standard Injection | LISR | 30 | 1 mm | 30° | 5 µm | |
| Standard Injection | LICR-TA35 | 35 | 1 mm | 30° | 5 µm | |
| Standard Injection | LICR-TA25 | 25 | 1 mm | 30° | 5 µm | |
| Standard Injection | LICR-TA20 | 20 | 1 mm | 30° | 5 µm | |
| Left Handed Injection | LISL | 30 | 1 mm | 30° | 5 µm | Left Handed |
| Short Spike Injection | LICR-SBS | 30 | 1 mm | 30° | 5 µm | Shorter Spike |
| Small ID Injection | LICR-ID4 | 30 | 1 mm | 24° | 4 µm | |
| Small ID Injection | LICR-ID4TA30 | 30 | 0.6 mm | 24° | 4 µm | |
| Small ID Injection | LICR-ID4TL0.6 | 35 | 0.6 mm | 24° | 4 µm | |
| Small ID Injection | LICR-OD5 | 30 | 1 mm | 20° | 4 µm | No Spike |

Holding Micropipettes

| Type | Reference no. | Available angles | Distance tip to bend | Inner diameter | Outer diameter |
|------------------|---------------|------------------|----------------------|----------------|----------------|
| Small Holding | LHC-OD75TA35 | 35 | 1 mm | 15 µm | 75 µm |
| Medium Holding | LHC-OD100TA35 | 35 | 1 mm | 15 µm | 100 µm |
| Medium Holding | LHC-OD100ID20 | 30 | 1 mm | 20 µm | 100 µm |
| Medium Holding | LHC-OD100TA20 | 20 | 1 mm | 20 µm | 100 µm |
| Standard Holding | LHC-ID15 | 30 | 1 mm | 15 µm | 120 µm |
| Standard Holding | LHS | 30 | 1 mm | 25 µm | 120 µm |
| Standard Holding | LHC-TA20 | 20 | 1 mm | 25 µm | 120 µm |
| Standard Holding | LHC-TA25 | 25 | 1 mm | 25 µm | 120 µm |
| Standard Holding | LHC-TA35 | 35 | 1 mm | 25 µm | 120 µm |
| Standard Holding | LHC-ID30 | 30 | 1 mm | 30 µm | 120 µm |
| Standard Holding | LHC-ID30TA35 | 35 | 1 mm | 30 µm | 120 µm |

Biopsy Micropipettes

| Type | Reference no. | Available angles | Distance tip to bend | Bevel | Inner diameter | Other |
|--------|---------------|------------------|----------------------|-------|----------------|---------|
| Biopsy | LBSR | 30 | 1 mm | 40° | 23 µm | Beveled |
| Biopsy | LBC-OD27BA90 | 30 | 1 mm | - | 21 µm | Flat |
| Biopsy | LBC-OD30BA90 | 30 | 1 mm | - | 23 µm | Flat |
| Biopsy | LBC-OD32TA35 | 35 | 1 mm | - | 25 µm | Flat |
| Biopsy | LBC-OD35BA90 | 30 | 1 mm | - | 27 µm | Flat |
| Biopsy | LBC-OD35TA35 | 35 | 1 mm | - | 27 µm | Flat |
| Biopsy | LBC-OD40BA90 | 30 | 1 mm | - | 31 µm | Flat |
| Biopsy | LBC-OD40TA35 | 35 | 1 mm | - | 31 µm | Flat |

Assisted Hatching Micropipettes

| Type | Reference no. | Available angles | Distance tip to bend | Inner diameter |
|---------------|---------------|------------------|----------------------|----------------|
| Zona Drilling | LDS | 30 | 1 mm | 9 µm |
| Zona Drilling | LDC-TA35 | 35 | 1 mm | 9 µm |

Zona Dissection Micropipettes

| Type | Reference no. | Available angles | Distance tip to bend |
|--------------|---------------|------------------|----------------------|
| Zona Cutting | LCS | 30 | 1 mm |
| Zona Cutting | LCC-TA35 | 35 | 1 mm |

The Pipette Company

Unit 13, 22 Ware Street,
Thebarton,
South Australia 5031
Tel: +61 8 8152 0266
www.pipetteco.com



CooperSurgical Companies