



# SYNGA FLEX

---

## Denudation and Manipulation System for Streamlined Routines in IVF Labs



## Synga Flex - reliable solution for your IVF lab

Comprehensive Range of Micropipettes for micromanipulation, denuding and handling:

- Close-tolerance and high quality medical grade plastic tips
- Flexible with smooth edges and consistent geometries - will not scratch Petri dishes
- Individually packed for sterility and supplied in multipacks for convenience and economy

# Denudation & Manipulation - Synga Flex

Synga Flex system offers comprehensive range of pipettes which will give you optimal control in micromanipulation. Pipettes of internal diameter 130 – 250 µm are suitable for denudation as well as cell manipulation. Strict quality control ensures the highest quality and consistency of all products. Finished product is sterilized by gamma irradiation and MEA tested.



Price-effective and user-friendly system. Colour coding of straws allows for straightforward differentiation of diameters. Pipettes are packaged in straws easily cut open by the Synga SG Cutter.

Order information: Synga Flex SG pipettes

Ref No.	Description	Inner	Package	Colour
SF-130-50	Denudation SG Pipette	Ø 130 µm	5 × 10 pipettes	red ●
SF-140-50	Denudation SG Pipette	Ø 140 µm	5 × 10 pipettes	green ●
SF-150-50	Denudation SG Pipette	Ø 150 µm	5 × 10 pipettes	white ○
SF-180-50	Manipulation SG Pipette	Ø 180 µm	5 × 10 pipettes	yellow ●
SF-250-50	Blastocyst SG Pipette	Ø 250 µm	5 × 10 pipettes	blue ●

**Products and innovations from experienced embryologists**  
We proudly serve life-focused laboratories.

